## Why you for having chosen VELPI!

VELPI has offered to professionals in the sector a range of sophisticated and reliable equipment, works according to ISO 9001, ISO 14001 and OHSAS 18001 Quality System Certification. Our instruments are built according to the International norms 10-1 and to the rules of the CE mark.

Product lines:

**Feed Line**  
DKL Digesters  
Recirculating water pump for fumes aspiration  
Circulator

Distillation units  
Solvent extractors  
Raw fiber extractors  
Fiber extractors  
Oxidation test reactor  
01 Elemental analyzer  
Consumables

**Environment Line**  
Thermoreactors  
Determination systems  
Preheated thermostats and incubators  
Floculators

Lead mixer  
Mineralization unit for trace heavy metals determination  
Turbidimeter  
Ion detector

**g Line**  
Magnetic stirrers  
Digital thermoregulator  
Magnetic stirrer  
Magnetic stirrers  
Magnetic stirrers  
Magnetic plates  
Magnetic mixers  
Magnetic generator

**Lab Solution**  
Circulating bath  
Vacuum pump  
Peristaltic pump

VELPI SCIENTIFICA srl  
Via Usmate (MB) Italy  
Via Usmate, 16  
Tel: +39 039 628811  
Email: [inse@velp.it](mailto:inse@velp.it)  
[www.velp.com](http://www.velp.com)

## Why you for having chosen VELPI!

VELPI offers to operators in the sector a wide range of sophisticated and reliable instruments, works according to the norms of the Certification del Sistema Qualità ISO 9001, ISO 14001 and OHSAS 18001. Our instruments are built according to the International norms 10-1 and to the rules of the CE mark.

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Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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
Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

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The heating plate type RC is designed for general use, for all those applications that require the heating of samples in specific containers positioned on the heating plate of the instrument.  
The temperature of the heating plate is controlled by an analogical bulb thermoregulator fitted directly into the heating plate. This makes it possible to obtain a constant control of the heating plate temperature from 50 to 370 °C.

**NOTE:** The container of the product used must be compatible with the temperature used. Normally, in these cases, Pyrex glass containers are used. The use of the heating plate with high temperatures may cause variations of the superficial colour but this doesn't change the characteristics of thermal, mechanical and chemical resistance.

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- RC Heating Plate
- Power supply cord
- Instruction manual

### 2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

### 2.2 Start-up

Rotate the temperature knob completely to the left. Then, set the temperature by turning the dedicated knob.

## 3. Operating controls

### TEMPERATURE KNOB

To select the temperature of the heating plate turn the corresponding knob (Heating).  
The temperature can be programmed between 50 and 370 °C.

### ON-OFF SWITCH

The on-off switch turns the unit on and off. If the switch is in the "OFF" position the unit is off; if the switch is in the "ON" position the unit is on.  
Always turn the unit off after use.

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.  
The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-inflammable non-aggressive detergent.

## 5. Technical data

Power supply	F20700174 230V/50-60Hz; F20710174 115V/50-60Hz
Dimensions WxHxD	165x115x280 mm
Weight	1.4 Kg
Overall power	600 W
Construction material	Epoxy painted aluminum structure
Diameter of the heating plate	155 mm
Programmable temperature range	50 – 370 °C
Type of temperature control	Analog
Overtemperature protection	Yes
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Level of electrical protection CEI EN60529	IP 42
Pollution degree CEI EN61010-1	2

1071 Hemispheric bowl for 250 ml flasks  
 1072 Hemispheric bowl for 500 ml flasks  
 1073 Hemispheric bowl for 1000 ml flasks  
 1069 Support rod

69 Fuse 5x20mm 5A (for code F20700174)  
 70 Fuse 5x20mm 8A (for code F20710174)  
 39 Foot 13Dx5H embedded  
 42 Thermostat knob 24D blue

		
A00001071	A00001072	A00001073
A00001069		

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